

## MACS Quant analyzer 16

	peak/window size (nm)	channel name	example fluorophore
405 Violet	450/50	V1	BV421
	525/50	V2	BV510
	579/34	V3	BV570
	615/20	V4	BV605
	667/30	V5	BV650
488 blue	525/50 nm	B1	FITC
	579/34 nm	B2	PE
	615/20	B3	PI
	725/40	B4	PerCP
	785/62	B5	Pe-Cy 5.5
	785/62	B6	Pe-Cy7
Red	667/30	R1	APC
	725/40	R2	AF700
	785/62	R3	APC-Cy7

For more information

<https://www.miltenyibiotec.com/US-en/products/macsqquant-analyzer-16.html?countryRedirected=1#gref>